



CY 2022 Real World Testing Plan for SecureEMR+

Executive Summary

This is the real world test plan for CY 2022 for SecureEMR+ certified EHR solution. It provides the real world test measurements and metrics that meet the intent and objectives of ONC's Condition of Certification and Maintenance of Certification requirement for real world testing (§ 170.405 Real world testing) to evaluate compliance with the certification criteria and interoperability of exchanging electronic health information (EHI) within the care and practice setting which it is targeted for use.

As ONC has stated in its rule, "The objective of real world testing is to verify the extent to which certified health IT deployed in operational production settings is demonstrating continued compliance to certification criteria and functioning with the intended use cases as part of the overall maintenance of a health IT's certification." We have worked toward this objective in designing our test plan and its subsequent real world testing measurements and metrics.

This document builds toward the final testing measurements and metrics we will use to evaluate our product interoperability within production settings. Within each measure, we document our testing methodology, the associated ONC criteria, justification for measurement, expected outcomes from the testing, care settings applied for this measure, and if applicable the number of clients to use in our real world testing approach.

We have included our timeline and milestones for completing the real world testing in CY 2022, and information about compliance with the Standards Version Advancement Process updates.

A table of contents is provided for quick access to any document section, including the testing measurements and metrics found at the end of this document. Our signed attestation of compliance with the real world testing requirements is on the following page.



Developer Attestation

This Real World Testing plan is complete with all required elements, including measures that address all certification criteria and care settings. All information in this plan is up to date and fully addresses the health IT developer's Real World Testing requirements.

Authorized Representative Name: Camilo Torres

Authorized Representative Email: ctorres@inmediata.com

Authorized Representative Phone: 787-783-3233

Authorized Representative Signature: *Camilo Torres*

DATE : 11/30/2021

SECUREEMR+

Executive Summary.....	1
Developer Attestation.....	2
General Information	4
Timeline and Milestones for Real World Testing CY 2022.....	5
Standards Version Advancement Process (SVAP) Updates	6
Real World Testing Measurements	7
Testing Methodologies	7
Number of Clients Sites	7
Care and Practice Settings Targeted.....	8
RWT Measure #1. Number of Transition of Care C-CDAs Successfully Sent	9
RWT Measure #2. C-CDAs Received, Processed, and Incorporated	11
RWT Measure #3. Number of NewRx and CancelRX Messages Successfully Sent	13
RWT Measure #4. eRx capabilities for controlled substances.....	15
RWT Measure #5. Number of Patient Batch Exports Run	16
RWT Measure #6. Number of Quality Measures Successfully Reported on to CMS.....	18
RWT Measure #7. Patient Portal Use.....	20
RWT Measure #8. Immunization Messages Successfully Sent to IIS/Immunization Registries..	22
RWT Measure #9. Number of Patient Immunization History Queries Sent	23
RWT Measure #10. Syndromic Surveillance Messages Successfully Sent to Public Registry	25
RWT Measure #11. How many different health care survey registries do you connect with? ..	26
RWT Measure #12. How many different applications/3rd party systems are using your API ...	27
capabilities?	
RWT Measure #13. Direct Messaging Capabilities	29



General Information

Plan Report ID Number: SecureEMR+_RWT_2022

Developer Name: Bizmatics Inc.

Product Name(s): SecureEMR+

Version Numbers(s): Denali 3.1

Certified Health IT Criteria: 315(b)(1)-(3), (b)(6); (c)(1)-(4); (e)(1); (f)(1)-(2), (7); (g)(7)-(9); (h)(1)

Product List (CHPL) ID(s) and Link(s):

- 15.04.04.2689.Prog.31.00.1.170929
- <https://chpl.healthit.gov/#/listing/8856>

Developer Real World Testing Page URL: <https://prognocis.com/macra/#rwt>

Timeline and Milestones for Real World Testing CY 2022

- 1Q-2022: Begin communication with clients to ask for their support and participation in real world testing. The goal is to have a sufficient number of clients committed for real world testing by the end of 1Q-2022.
- 2Q-3Q 2022. During the 2nd and 3rd quarter of CY 2022, the real world testing with clients will be scheduled and performed. It is expected that a preparatory call will be done with clients to prepare them for testing activities. Results will be documented in the test results section of the test methods and ultimately used to build the test report. If any non-compliances are observed, we will notify the ONC-ACB of the findings and make the necessary changes required.
- 4Q-2022. During the last quarter of the year, the CY 2023 real world test plan will be completed according to ONC and ONC-ACB requirements and expectations. Test plan will be prepared for submission before the end of the year.
- 1Q-2023. Submit RWT Test Report to ONC-ACB.

Standards Version Advancement Process (SVAP) Updates

For CY 2022, we are not planning to make any version updates on approved standards through the SVAP process. We plan on implementing USCDI v1 in our C-CDAs and API support during CY 2022, but we have not finalized an exact date for rollout.

Standard (and version)	N/A
Updated certification criteria and associated product	N/A
Health IT Module CHPL ID	N/A
Method used for standard update	N/A
Date of ONC-ACB notification	N/A
Date of customer notification (SVAP only)	N/A
Conformance measure	N/A
USCDI-updated certification criteria (and USCDI version)	N/A

Real World Testing Measurements

The measurements for our real world testing plan are described below. Each measurement contains:

- Associated ONC criteria
- Testing Methodology used
- Description of the measurement/metric
- Justification for the measurement/metric
- Expected outcomes in testing for the measurement/metric
- Number of client sites to use in testing (if applicable)
- Care settings which are targeted with the measurement/metric

In each measurement evaluate, we elaborate specifically on our justification for choosing this measure and the expected outcomes. All measurements were chosen to best evaluate compliance with the certification criteria and interoperability of exchanging electronic health information (EHI) within the certified EHR.

Testing Methodologies

For each measurement, a testing methodology is used. For our test plan, we use the following methodologies.

Reporting/Logging: This methodology uses the logging or reporting capabilities of the EHR to examine functionality performed in the system. A typical example of this is the measure reporting done for the automate measure calculation required in 315(g)(2), but it can also be aspects of the audit log or customized reports from the EHR. This methodology often provides historical measurement reports which can be accessed at different times of the year and evaluate interoperability of EHR functionality, and it can serve as a benchmark for evaluating real world testing over multiple time intervals.

Compliance and/or Tool: This methodology uses inspection to evaluate if EHR is compliant to the ONC criteria requirements. It can be done through 1-v-1 inspection testing or utilize various tools to measure or evaluate compliance and interoperability. If an EHR Module capabilities is not widely used in production by current users, compliance inspection can provide assurance criteria is working as previously certified.

Survey/Self-Test: This methodology evaluates interoperability and compliance of EHR Module capabilities through feedback from users. ONC has recognized that self-testing can be a viable method for evaluation and compliance, and this methodology can provide insight into how clinicians employ and use a feature which reveals actual value and impact of interoperability of the EHR Module.

Number of Clients Sites

SECUREEMR+

Within each measure, we note the minimum number of clients or client sites we plan to use for this measure evaluation. The numbers vary depending on the methodology as well as overall use of the associated EHR Module criteria by our users. For criteria that are not widely used by our customer base, we may test the respective measure in our own production-sandbox environment given lack of customer experience with the criteria functionality.

Care and Practice Settings Targeted

Our EHR targets different care settings, including family practice, internal medicine, general practice, pediatrics, and some other settings. While we do have different physician types using our EHR, our workflow design is based on general ambulatory encounters. In each measure, we do also address the care settings targeted and note any necessary adjustment or specific factor to consider with this specific measure.

RWT Measure #1. Number of Transition of Care C-CDAs Successfully Sent

Associated Criteria: 315(b)(1)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking and counting how many C-CDAs are created and successfully sent from the EHR Module to a 3rd party via Direct messaging during a transition of care event over the course of a given interval.

Measurement Justification

This measure will provide a numeric value to indicate both the how often this interoperability feature is being used as well as its compliance to the requirement. An increment to this measure indicates that the EHR can create a C-CDA patient summary record, including ability to record all clinical data elements, and by sending the C-CDA patient summary record, the EHR demonstrates successful interoperability of an exchanged patient record with a 3rd party. This measurement shows support for Direct Edge protocol in connecting to a HISP for successful transmission.

Measurement Expected Outcome

The measurement will produce numeric results over a given interval. We will utilize various reports and audit logs, including Automated Measure (315.g.2) reports, to determine our measure count. During the year, we will examine the log information for a minimum period of three (3) months to determine an appropriate sample of this measurement.

A successful measure increment indicates compliance to the underlying ONC criteria. It will show that the EHR can create the C-CDA patient summary record, including record required clinical data elements. In sending the C-CDA patient summary record, the EHR will demonstrate ability to confirm successful interoperability of an exchanged patient record with a 3rd party, including support for Direct Edge protocol in connecting to a HISP. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

We will use the measure count to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

SECUREEMR+

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s). This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

RWT Measure #2. C-CDAs Received, Processed, and Incorporated Associated Criteria: 315(b)(1) and 315(b)(2)

Testing Methodology: Reporting/Logging and Compliance

Measurement Description

This use case is associated with two different criteria and tracking two different metrics. It is tracking and counting how many C-CDAs are successfully received and/or incorporated upon receipt from a 3rd party via Direct messaging during a transition of care event over the course of a given interval. It is also doing a compliance evaluation of the inbound C-CDA error detection capability.

Measurement Justification

For the first measure metric, we will provide a numeric value to indicate both the how often the EHR can receive a C-CDA patient summary record and then incorporating problems, medications, and medication allergies from the C-CDA into the patient record. This measure also shows support for Direct Edge protocol in connecting to a HISP for successful transmission.

The second measure is evaluation the C-CDA error detection capability. C-CDA error detection provides assurance to the user of the validity of received or imported in C-CDAs which is both a certification requirement and supports interoperability within production setting.

Measurement Expected Outcome

This use case has two measures.

Measure #1: C-CDA Incorporation. The first measurement will produce numeric results over a given interval. We will utilize various reports and audit logs, including Automated Measure (315.g.2) reports, to determine our measure count. During the year, we will examine the log information for a minimum period of three (3) months to determine an appropriate sample of this measurement.

A successful measure increment indicates compliance to the underlying ONC criteria. It will show that the EHR can create the EHR can receive a C-CDA patient summary record. In incorporating the C-CDA patient summary record, the EHR will demonstrate successful interoperability of problems, medications, and medication allergies of patient record with a 3rd party, including support for Direct Edge protocol in connecting to a HISP. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

SECUREEMR+

We will use the measure count to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

Measure #2: C-CDA Error Detection. SecureEMR+ has an internal tool for tracking the errors from C-CDA received. The user will import in, either through upload or inbound messages, C-CDA's with different known errors. The user will use the EHR functions to parse the C-CDA document and perform errors detection which will be reviewed by the user. We will confirm the process and steps done by the user meet the criteria requirements of the EHR Module and works as expected in production-type environment.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s). This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

RWT Measure #3. Number of NewRx and CancelRX Messages Successfully Sent

Associated Criteria: 315(b)(3)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking and counting how many NewRx and CancelRx electronic prescriptions were created and successfully sent from the EHR Module to a pharmacy destination over the course of a given interval.

Measurement Justification

This measure will provide a numeric value to indicate both the how often this interoperability feature is being used as well as its compliance to the requirement. An increment to this measure indicates that the EHR can create a NewRx SCRIPT electronic prescription message as well as the CancelRx message and transmit it to a pharmacy, typically via the Surescripts Network.

Measurement Expected Outcome

This use case will have two different measurement metrics, and both measurements will produce numeric results over a given interval. We will utilize various reports and audit logs, including Automated Measure (315.g.2) reports, to determine our measure count. During the year, we will examine the log information for a minimum period of three (3) months to determine an appropriate sample of this measurement.

Measure #1 will show that the EHR can create the NewRx message and send over a production network, like the Surescripts Network, to a pharmacy. Measure #2 will count the number of CancelRx message sent over the same time period.

Successfully completing these measure use cases also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

We will use the measure counts to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

SECUREEMR+

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s). This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

RWT Measure #4. eRx capabilities for controlled substances

Associated Criteria: 315(b)(3)

Testing Methodology: Reporting

Measurement Description

This is a measure to determine how often providers are electronically prescribing controlled medications from the EHR.

Measurement Justification

This measure will evaluate how often provider electronically prescribed controlled substance medications from their EHR.

Controlled substances are not an explicit requirement of ONC certification, but the electronic prescribing features should support controlled substances if the provider and support health IT system have the other capabilities enabled for electronic prescribing of controlled substances.

Measurement Expected Outcome

We will access our eRx system to run detailed reports on providers who use EPCS. We will count the controlled substances successfully electronically prescribed. If we observe any errors, we will investigate and address accordingly.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. During the year, we will examine the log information for all providers for a minimum period of one (1) month to determine how many controlled substances our system electronically prescribes for our provider community.

RWT Measure #5. Number of Patient Batch Exports Run

Associated Criteria: 315(b)(6)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking and counting how many batch exports of C-CDAs were successfully performed by the EHR Module over the course of a given interval.

Measurement Justification

This measure will provide a numeric value to indicate both the how often this interoperability feature is being used as well as its compliance to the requirement. An increment to this measure indicates that the EHR can create a batch export of multiple C-CDA patient summary records.

Measurement Expected Outcome

The measurement will produce numeric results over a given interval. We will utilize various reports and audit logs to determine our measure count. During the year, we will examine the log information for a minimum period of three (3) months to determine an appropriate sample of this measurement.

A successful measure increment indicates compliance to the underlying ONC criteria. It will show that the EHR can create a batch export of multiple C-CDA patient summary records, which can be used in means of health IT interoperability. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

We will use the measure count to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s). This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

SECUREEMR+

RWT Measure #6. Number of Quality Measures Successfully Reported on to CMS

Associated Criteria: 315(c)(1)-(c)(4)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking and counting how many eCQM quality measures were successfully reported on by the EHR Module to CMS during their submission period for MIPS Quality reporting.

Measurement Justification

This measure will provide a count and list of electronic clinical quality measures (eCQMs) which are calculated and submitted to CMS for a given program, like MIPS. Clinical quality measures are only used for the respective CMS programs and any production measures should utilize submission to CMS. Because CQM criteria, 315(c)(1)-(c)(4), all work collectively together in the eCQM functionality of the EHR Module, this measurement is used for all three.

Measurement Expected Outcome

The measurement will a count and list of eCQMs submitted to CMS over a given interval. We will utilize our population health tool to determine our measure count.

A successful measure submission indicates compliance to the underlying ONC criteria. It will show that the EHR can do calculations on the eCQM and that they are accepted by CMS. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

We will use the measure result to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s). This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

SECUREEMR+

RWT Measure #7. Patient Portal Use

Associated Criteria: 315(e)(1)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking and counting how many patients are successfully logged into and accessed their patient portal account as well as email transmissions from the portal over the course of a given interval.

Measurement Justification

This use case has three metrics which measure real world interoperability actions of patients and their secure portals. The first measure will provide a numeric value to indicate both the how often patients log into their patient portal to view, download, or transmit their health data. The second measure counts the number of email transmission will be done from the portal, and the third measure tracks how many different destinations transmission are sent.

These activities show how patients commonly utilize a patient portal as well as the breadths of its use with other health care entities.

Measurement Expected Outcome

The measurements will produce numeric results over a given time interval of a minimum of three (3) months. We will utilize various reports and audit logs, including Automated Measure (315.g.2) reports, to determine our measure count.

For all three measures, a successful increment indicates compliance to the underlying ONC criteria. It will show that patients can log into their patient portal to access their patient data and transmitting their health data to a 3rd party. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

We will use the measure counts to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s).

SECUREEMR+

This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

RWT Measure #8. Immunization Messages Successfully Sent to IIS/Immunization Registries

Associated Criteria: 315(f)(1)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking how immunization messages are created and successfully sent from the EHR Module to an IIS/immunization registry over the course of a given interval.

Measurement Justification

This measure will provide a value to indicate if bi-direction exchange is enabled and working between the EHR and the IIS/immunization registry. This measure indicates that the EHR can exchange an immunization message, including ability to record all clinical data elements, with an IIS/immunization registry.

Measurement Expected Outcome

The measurement will record if the immunization interface is active and working between the EHR and public health registry. Our IIS interface log/report indicates status of the interface, and during the year, we will examine the log information for a minimum period of three (3) months to determine an appropriate sample of this measurement.

A successful measurement status indicates compliance to the underlying ONC criteria. It will show that the EHR is creating and exchange the HL7 immunization record, including ability to record the required clinical data elements. In sending the immunization message, the EHR will demonstrate ability to confirm successful interoperability of patient's immunization data to an IIS/immunization registry. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s). This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

RWT Measure #9. Number of Patient Immunization History Queries Sent

Associated Criteria: 315(f)(1)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking and counting how many patient immunization history queries are created and successfully sent from the EHR Module to an IIS/immunization registry over the course of a given interval.

Measurement Justification

This measure will provide a numeric value to indicate both the how often this interoperability feature is being used as well as its compliance to the requirement. An increment to this measure indicates that the EHR can create a patient immunization history query message, and by sending the message, the EHR demonstrates successful interoperability with an IIS/immunization registry.

Measurement Expected Outcome

The measurement will produce numeric results over a given interval. We will utilize various reports and IIS transaction logs to determine our measure count. During the year, we will examine the log information for a minimum period of three (3) months to determine an appropriate sample of this measurement.

A successful measure increment indicates compliance to the underlying ONC criteria. It will show that the EHR can create the HL7 immunization history query of a patient. In sending the immunization message, the EHR will demonstrate ability to confirm successful interoperability with an IIS/immunization registry. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

We will use the measure count to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this

SECUREEMR+

measure to be applicable to all settings. We will test a minimum of five (5) client practice(s). This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.

RWT Measure #10. Syndromic Surveillance Messages Successfully Sent to Public Registry

Associated Criteria: 315(f)(2)

Testing Methodology: Reporting/Logging

Measurement Description

This measure is tracking how many syndromic surveillance messages are created and successfully sent from the EHR Module to a syndromic registry over the course of a given interval.

Measurement Justification

This measure will provide a value to indicate if the exchange is enabled and working between the EHR and the syndromic surveillance registry. This measure indicates that the EHR can exchange a syndromic message, including ability to record all clinical data elements, with the public registry.

Measurement Expected Outcome

The measurement will record if the syndromic surveillance interface is active and working between the EHR and public health registry. Our logs indicate the status of the interface, and during the year, we will examine this log information for a minimum period of three (3) months to determine an appropriate sample of this measurement.

A successful measure status indicates compliance to the underlying ONC criteria. It will show that the EHR can create the HL7 syndromic surveillance message, including ability to record the required clinical data elements. In sending the syndromic surveillance message, the EHR will demonstrate ability to confirm successful interoperability of patient's syndromic data to public health registry. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will survey a sample of our client community to obtain our survey results.

RWT Measure #11. How many different health care survey registries do you connect with?

Associated Criteria: 315(f)(7)

Testing Methodology: Survey/Self-Test

Measurement Description

This is a survey measure to determine the number of health care survey public health registries you use.

Measurement Justification

This measure will survey users to determine real world interoperability and usability, specifically many different public health care survey registries are used by the provider.

A survey or self-testing can often provide more information on the impact and value of an interoperability element than a standard software test evaluation. This survey measure will the number and names of health care survey registries which are integrated with the EHR.

Measurement Expected Outcome

The user will be asked the survey question and given the survey answer choices below:

- Numeric answer to the question, and if willing, the names of the other systems.

The answer will provide insight into how clinicians view both the use and value of this interoperability feature. For example, response may show that additional training is needed to better utilize the feature or that it is not currently utilized as currently designed. It will provide a benchmark for evaluate future surveys as well as to share insight into any new development for improvements or enhancements of the health IT system.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will survey a sample of our client community to obtain our survey results.

RWT Measure #12. How many different applications/3rd party systems are using your API capabilities?

Associated Criteria: 315(g)(7)-(g)(9)

Testing Methodology: Survey/Self-Test

Measurement Description

This is a survey measure to determine how many different systems or applications are connecting to your EHR via the API.

Measurement Justification

This measure will survey users to determine real world interoperability and usability, specifically many 3rd party systems or applications are integrated and using the EHR's API interface.

We do not believe many of our clients are using API capabilities and therefore do not believe a measurement metric will be sufficient. Instead, we believe best means to do real world testing in this initial year of the program is to survey users to determine API usage and then factor that information into future real world testing efforts.

A survey or self-testing can often provide more information on the impact and value of an interoperability element than a standard software test evaluation. API capabilities are an important component of the modern health IT system, and utilization of API resources will help improve patient care and care coordination.

Measurement Expected Outcome

The user will be asked the survey question and given the survey answer choices below:

- Numeric answer to the question, and if willing, the names of the other systems.

The answer will provide insight into how clinicians view both the use and value of this interoperability feature. For example, response may show that additional training is needed to better utilize the feature or that it is not currently utilized as currently designed. It will provide a benchmark for evaluate future surveys as well as to share insight into any new development for improvements or enhancements of the health IT system.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will survey a sample of our client community to obtain our survey results.

SECUREEMR+

RWT Measure #13. Direct Messaging Capabilities

Associated Criteria: 315(h)(1)

Testing Methodology: Reporting/Logging

Measurement Description

This use case is tracking various measurements associated with Direct messaging successfully sent and received with the EHR Module and a 3rd party over the course of a given interval.

Measurement Justification

This use case has three related measure metrics. The first is number of Direct messages sent. The second is the number of different 3rd party destinations the Direct messages were sent to, and the third measure is the number of Direct messages received by the EHR from a 3rd party. These measures will reveal the volume and depth of Direct messaging done in real world.

Measurement Expected Outcome

The measurements will produce numeric results over a given interval. We will utilize various reports and audit logs, including Automated Measure (315.g.2) reports, to determine our measure count.

For each of the three measures, a successful increment indicates compliance to the underlying ONC criteria. It will show that the EHR can be authenticated with DirectTrust, create a Direct message, and demonstrate interoperability of an exchanged message with a 3rd party. Successfully completing this measure also implies users have a general understanding of the EHR functional operations for this EHR Module and an overall support for the user experience while not completing this measure may indicate lack of understanding or possibly lack of use or need for this functionality. For the second measure, a higher number indicates the interoperability feature is utilized across a wide range of different partners while a smaller number indicates a more focused distribution.

We will use the measure count to establish a historic baseline of expected interoperability use so it can be used in subsequent real world testing efforts.

Care Settings and Number of Clients Site to Test

Our EHR is targeted to a variety of care settings including family practice, internal medicine, general practice, pediatrics, and some other settings in ambulatory care. We have designed this measure to be applicable to all settings. We will test a minimum of five (5) client practice(s).

SECUREEMR+

This number covers a sufficient percentage of existing practices to provide a viable sample of users of the certified EHRs.